

Strap-on Intramuscular Injection Trainer

For the practice between students, simulated patients and mannikin, facilitating their understanding of the intramuscular injection process and communication.



FEATURES

- 1 | The water tank is located inside this trainer.
- 2 | This trainer is easier to drain without the waste water bag.
- 3 | Realistic tactile feeling of skin and acromion allows palpating the anatomical landmark.

SKILLS

- | Intramuscular Injection training
- | COVID-19 vaccination training

SPECIFICATIONS

Size: W10.5×D26×H5.5cm
W4.2 x D10.3 x H2.2 in.

SET INCLUDES

- 1 main unit
(frame, skin, soft tissue, muscle,
acromion, drainage reservoir, belts)
- 1 instruction manual

DESCRIPTIONS

MATERIALS

Soft resin
Latex free

RECOMMENDED DEVICES

IM needle 23G

REPLACEMENT PARTS

- 11297-220 replacement skin
- 11297-230 replacement soft tissue part
- 11297-240 replacement muscular part