

# Airway Suction Trainer

# Airway suction and tracheostomy management simulator

Life-like anatomical and tactile features to further understanding temporary airway suctioning.





Nasal suction



Oral suction

## **FEATURES**

- 1 | Facilitate understanding of the way to insert catheter.
- 2 | Transparent cross-sectional window to monitor the catheter
- 3 | Realistically viscous simulated sputum is included
- 4 | Can be easily dismantled for cleaning
- 5 | Able to open the mouth by hand

## **SKILLS**

- Oral suctioning
- | Nasal suctioning
- | Endotracheal suctioning

# LANDMARKS

- oral cavity
- tongue
- oral and nasal pharynx
- | epiglottis
- | trachea
- | esophagus



## 11229-060 Oral Cover for training using liquid

Oral cover to store simulated sputum can be purchased as optional parts. This oral cover is to maintain the sputum inside the oral cavity.



## **DESCRIPTIONS**

#### SPECIFICATIONS

Size: 45x23x15cm / 17.7x9x5.9 inch Weight: 2kg/4.4lbs

#### SET INCLUDES

- 1 head & neck model with airway anatomy bottle of lubricant 2 catheters
- bottle of simulated sputum
- storage case lubricant instruction manual

1 syringe

#### **MATERIALS**

Soft resin Latex free

### REPLACEMENT PARTS

OPTIONAL PARTS 11229-020 1 replacement skin with trachea 11229-060 oral cover for training using liquid

- 11229-070 1 simulated sputum
- 11229-040 2 sample catheters 11229-050 1 lubricant





