



MW64

Product Code: 11454-999

Vaccination Training Set (Set of 5)

MW 64 Vaccination Training Set will undeniably help in training the thousands of healthcare workers who are involved in COVID-19 vaccination because it is simple, easy to use, and low maintenance.



Acromion
(Landmark)

Depth: 1 1/2"



The hole for
draining the water

FEATURES

- | The pad simulates deltoid muscle of arm.
- | The pad's depth is 1 1/2" (38.1mm)
- | This vaccination training set can be used together with the Kyoto Kagaku patient care mannequins and students.

SET INCLUDES

- 5 Vaccination Training Pads/
- 5 Pad holders/ 5 Strap-on belts/
- 1 Instruction manual

REPLACEMENT PART

- Pad for Vaccination Training Set (Set of 5)
- 11454-010

 **KYOTO KAGAKU co.,LTD**

 <http://www.kyotokagaku.com>
 rw-kyoto@kyotokagaku.co.jp

**Worldwide Inquires and Ordering
Head Office and Factories**

15 Kitanechoya-cho, Fushimi-ku,
Kyoto, 612-8388, Japan

 +81-75-605-2510  +81-75-605-2519

Hong Kong Representative Office

Level 60, One Island East,
18 Westlands Road, Island East, Hong Kong

 +852-3978-1365  +852-3978-1199

European Representative Office

De-Saint-Exupery-Str.10,
60549 Frankfurt, Germany


 +49-69-5060-28160  +49-69-5060-28170

 All American regions
KYOTOKAGAKU AMERICA INC.

 <http://www.kkamerica-inc.com/>
* United States only
 info@kkamerica-inc.com

3109 Lomita Boulevard, Torrance,
CA 90505-5108, USA

 +1-310-325-8860

 +1-310-325-8867