

Casualty Care Rescue Randy product profile

Product Name	Casualty Care Rescue Randy - powered by Strategic Operations Hyper-Realistic® technology
Addresses Training Needs/ Curricula	TCCC - Tactical Combat Casualty Care for all service members TECC - Tactical Emergency Casualty Care for all EMS/EMT/LEO professionals
Configuration	<ul style="list-style-type: none"> • Full-body manikin with articulated joints • Weight-distributed to mimic the human anatomy • Holds 3-4 liters of blood • Pumping system with regulator simulating human blood flow pressure (2-3) psi • Manual foot operated blood pumping system allows 3 wounds to bleed simultaneously • Optional Electric Deluxe blood pumping system allowing 4 wounds to bleed simultaneously • Comes with blood powder (makes 5 gallons) • Three (3) repairable neck skins • Ten (10) repairable tracheas • One (1) repair kit
Dimensions	Male at 72" tall x 25" wide x 12" deep (165.1 cm x 63.5 cm x 30.48 cm)
Product Ship Weight	122 lbs. (2 boxes)
Product Code/SKU	149-6075 - Casualty Care Rescue Randy 149-6095 - Casualty Care Rescue Randy with Basic Pumping System (BPS)
Consumables	<ul style="list-style-type: none"> • 149-6079 Trachea, Randy - (sold in packs of 10 qty) • 149-6080 Repairable neck skin Randy - (sold in packs of 3 qty) • 149-6081 Silicone Repair Kit - (1 kit)
Additional Options	<ul style="list-style-type: none"> • IV and Suture Sleeve / - # IVSS • IV Veins (package of 10) / - # IVV



Nasco Healthcare Inc.
16 Simulaid Drive
Saugerties, NY, 12477 USA
Phone: 1-833-NASCOHC (627-2642)
Email: info@nascohealthcare.com
www.nascohealthcare.com



Ready to treat. Ready to act.

Introducing Casualty Care Rescue Randy

Full-body manikin that is uniquely designed for enhanced realistic training on the procedures that treat the 3 most preventable causes of death:

1. Massive Bleeding
2. Tension Pneumothorax
3. Airway Obstruction

- Compliant with TCCC/TECC scenarios
- Rugged, reliable, and hyper-realistic
- 100% hand-made in the USA

Available Blood Pumping System #BPS2-B


Allows for all four CCRR wounds to bleed simultaneously and includes additional wounds

Included:

- Backpack with pump and blood reservoir (1)
- Wireless remote key fob (2)
- Rechargeable batteries (2)
- Battery charger (1)
- Upper extremity bleeding module (1)
- Lower extremity bleeding module (1)



Casualty Care Rescue Randy:
Life-like training solutions for
real-life situations.



This manikin is intended to be used as an instructional aide in the teaching and demonstrating of live rescue training. Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.