

E-SCOPE®

HEARING IMPAIRED



Used with hearing aid or hearing impairment



One-handed control over volume adjustments



Specialist adult diaphragm



Automatic shut off after 1.5 to 2 minutes



Audio output jack



Battery operated (single AAA battery)




Supplied with accessory pack



The **E-SCOPE Hearing impaired Model** is a range of auditory accessories, allowing ability to choose the appropriate specific needs. This model is ideal for the **Clinical Model** — amplifying sound from an acoustic scope. Separate headset offered to ensure we find the right fit with hearing loss.

E-SCOPE Catalog # 718-7700 (Clinical model)
718-7710 (Hearing impaired model)
Technical Support: 011 792 2190
Warranty: One year



Survival Technology

Your Link in the
Chain of Survival

www.survivaltechnology.com

6 Susan Street, Strydom Park Randburg

info@survivaltechnology.com

Tel: +27 11 792 2190

Fax: +27 11 793 4234

PO Box 416, Ferndale, 2160